7.20-8.27

1.Constructed the plasmid gYB2a-pobRWT-mCherry-codA-cmr The person primarily responsible for the work: Qianwen Jin Construction time:

1. Obtain CD fragment by PCR

The target fragment is located on the pUAM-DE-CD plasmid, and the amplified target fragment is obtained by PCR. The PCR system and procedure are showed below.

PCR system	(50µL)	
pUAM-DE-CD	10ng	
CD-Gibson-0425 F	2μL	
CD-Gibson-0425 R	2µL	
2 x Mix	25µL	
DDW	to 50µL	
Table 1		

The PCR products are detected by agarose gel electrophoresis, and the correct target fragment is 1332bp. We obtained the correct target fragment, and the sample is purified.

2. Obtain gYb2a-PpobA*2-mCherry-Cmr fragment by PCR

The target fragment is located on the gYb2a-PpobA*2-mCherry-SacB-Cmr plasmid, and the amplified target fragment is obtained by PCR. The PCR system and procedure are showed below.

PCR system (50µL)		
gYb2a-PpobA*2-mCherry-		
SacB-Cmr	10ng	
Mc-Gibson-F	2µL	
Mc-Gibson-R	2µL	
2 x Mix	25µL	
DDW	to 50µL	
Table 2		

The correct target fragment is 4572bp. We obtained the correct target fragment, and the sample is purified.

3. Gibson connection

The CD fragment and gYb2a-PpobA*2-mCherry-Cmr fragment are connected by Gibson connection method, and the connection system is as follows.

Connection system (10µL)	
CD	2.4µL
gYb2a-PpobA*2-mCherry-	1μL
XX-Cmr	
5 x Cell Buffer	2μL
Exnase II	1µL
DDW	3.6µL
Table 3	

4. Colony PCR

-

After the petri dish is incubated at 37°C for 12 hours, 5 colonies were selected on the plate. The colony PCR system and procedure were as follows.

PCR system (10µL)		
gYb2a-PpobA*2-mCherry-	1µL	
CD-Cmr		
CD-Gibson-0420-F	0.4µL	
CD-Gibson-0420-R	0.4µL	
2 x Mix	5μL	
DDW	3.2µL	

Table 4

The PCR products were detected by agarose gel electrophoresis, and the results were as follows



Figure 1

The correct target fragment is about 1300bp, and the length of the colony PCR sample in lanes 3,4 and 5 is inferred from the gel electrophoresis image is correct. The above 3 strains were expanded and the plasmids were put forward.

5. Enzyme digestion verification

Use the enzyme KpnI to cut the plasmid gYb2a-PpobA*2-mCherry-CD-Cmr at the same time to

verify whether the plasmid is constructed correctly. The following is the system of digestion.

2 - goodon of soom (1	• • • • • •
gYb2a-PpobA*2-mCherry-	100ng
CD-Cmr	
KpnI	0.2µL
custsmart	1µL
DDW	to 10µL

Digestion system (10µL)

Table 5

Digested gYb2a-PpobA*2-mCherry-CD-Cmr is verified by electrophoresis which is showed below.



It is inferred from the gel electrophoresis that the length of the sample is as expected. The construction is preliminarily correct. We sequenced the constructed plasmid. The result shows that the sequence is the same as expected, indicating that our target plasmid has been constructed successfully.

8.28-9.3

2. Validation of PobRWT

The person primarily responsible for the work: Xiaoya Wei Experimental time:

PobR^{WT} exhibited a high degree of ligand specificity with a narrow dynamic range.

Original wild type colonies was induced with 10-6,10-5,10-4,10-3,10-2,0.1,1,10g/L 4HB, and after 12 incubations, the optical density at 600 nm (OD600) and red fluorescence were measured (552 nm as excitation wavelength and 600 nm as emission wavelength), and the fluorescence response curve was plotted.

Experimental group (three groups in parallel):

200µL M9 medium + 1mg/L Amp + 4HB at each concentration+ bacteria





9.5-9.7

3. Validating the feasibility of a dual selection system

The person primarily responsible for the work: Jianing Li Experimental time:

In our design, with the presence of 4HB, the *E. coli* BW $\Delta codA$ harboring the PobR^{WT} biosensor was resistant to chloramphenicol, but sensitive to 5-FC.

Our verification methods: Cultivation of bacteria harboring the PobR^{WT} biosensor in medium containing different concentrations of chloramphenicol and supplied with 0.5 g/L 4HB.

Cultivation of bacteria in medium containing either 0.5 g/L 4HB or 5 mg/L 5-FC alone, or both. The OD_{600} was measured after 12h of cultivation.



9.8-9.10

4. Test the capacity of our mutant library

We used the dilution coating plate method to determine the total library volume. We diluted 50ul of the bacterial solution to 5000ul (1:100), from which 50ul of the solution was aspirated and applied to a plate with a final concentration of 0.1mg/L Amp, and the total amount of the mutant library was estimated based on the number of colonies versus the dilution multiple.



The selection capacity for each compound was more than 900,000 clones (with at least four plates about 450,000 CFU/mL.)

9.5-9.10

5. Dual selection significantly reduces the response of the sensor to natural substrates As these mutants underwent two rounds of negative selections, when tested with increasing 4HB concentration, the obtained PobR mutants exhibited significantly reduced basal expression and robustly attenuated 4HB responsiveness compared to that of PobR^{WT}.



9.11-9.29

6.Selecting the PobR mutation library

In positive selection experiments of PobR mutants using different ligands, we obtained a biosensor that the fluorescence intensity was at least 1.5-fold higher than that of the negative control.

Take 500 μ L of the second-round screening bacterial solution and spread it on the solid medium supplemented with 0.1 mg/L of ampicillin, 0.09 g/L of chloramphenicol and 0.5 g/L of HPP /0.3 g/L of PPA /0.3 g/L of 2-PE /0.3 g/L of HMA /0.5 g/L of MA /0.3 g/L of PAld /0.3 g/L of p-Coumaric to isolate single colonies. Pick the red

colonies from the selection medium and culture them in liquid LB containing ampicillin for 8-10 hours.

After culturing for 8-10 hours, take 2 microliters of bacterial broth to inoculate 200 microliters containing ampicillin, 0.5 g/L of HPP, and ampicillin only. in M9 medium. After 12 hours, measure the optical density at 600 nm (OD600) and red fluorescence (552 nm as excitation wavelength and 600 nm as emission wavelength).

Induction group: 200 μL M9 medium + 0.1mg/L Amp + 0.5 g/L of HPP /0.3 g/L of PPA /0.3 g/L of 2-PE /0.3 g/L of HMA /0.5 g/L of MA /0.3 g/L of PA1d /0.3 g/L + bacterial broth.

Uninduced group: 200 μL M9 medium + 0.1mg/L Amp + bacterial broth.







Among these PobR mutants, the best responsive sensor was the PobR mutant clone 0914-A8-1 with a point mutation of W177R, which showed 5.1-fold increase of mCherry fluorescence in response to HPP.

7.Probing the specificity of PobR^{W177R} to substrates and the range of response

We assessed the specificity of PobRW177R in detecting different aromatic compounds, it was induced with uninduced,0.5g/L of HPP,4HB,HMA,PPA,2-PE,MA,PALD,L-Phe and p-coumaric acid,and after 12 incubations, the optical density at 600 nm (OD600) and red fluorescence were measured (552 nm as excitation wavelength and 600 nm as emission wavelength), and observed that 4HB could still trigger its activity in driving mCherry expression.

