

9.pRBI2s-single-T7-STS Xinyao Yuan

9.1 Enzyme digestion

Use two enzymes Xho I to cut the plasmid pRBI2s-single-T7 The following is the system of digestion system.

Enzyme digestion system (50μL)	
pRBI2s-single-T7	10ng
Xho I	1μL
Cutsmart	5μL
DDW	To 50μL

Table 1

The Enzyme digestion products are detected by agarose gel electrophoresis, and the correct target fragment is 3500bp. We obtained the correct target fragment, and the sample is purified.

9.2 Obtain STS by PCR

The target fragment is located on the STS plasmid, and the amplified target fragment is obtained by PCR. The PCR system and procedure are showed below.

PCR system (50μL)	
STS plasmid	10ng
STS-G-F	2μL
STS-G-R	2μL
2 x Mix	25μL
DDW	To 50μL

Table 2

The PCR products are detected by agarose gel electrophoresis, and the correct target fragment is 1179bp. We obtained the correct target fragment, and the sample is purified.

9.3 Gibson connection

The STS fragments and pRBI2s-single-T7 are connected by Gibson connection method, and the connection system is as follows.

Connection system (20μL)	
pRBI2S-single-T7	1.2μL
STS	1.4μL
5×CE Buffer	2μL
ExhaseII	1μL

DDW	4.4 μ L
-----	-------------

Table 3

9.4 Colony PCR

After the petri dish is incubated at 37°C for 12 hours, 11 colonies were selected on the plate. The colony PCR system and procedure were as follows.

PCR system (10 μ L)	
pRBI2s-single-T7-STS	1 μ L
STS-JUNP-F	0.2 μ L
RSF-JUNP-R	0.2 μ L
2 x Mix	5 μ L
DDW	to 10 μ L

Table 4

The PCR products were detected by agarose gel electrophoresis, and the results were as follows.




Figure 1

The length of the colony PCR sample in lanes 3、5、7、8、9、13 are inferred from the gel electrophoresis image are correct. The 7 and 8 strains were expanded and the plasmids were put forward.

9.5 Enzyme digestion verification

Use Xho I to cut the plasmid pRBI2s-single-T7-STS to verify whether the plasmid is constructed correctly. The following is the system of digestion.

Digestion system (10 μ L)	
pRBI2s-single-T7-STS	100ng
Xho I	0.2 μ L
cutsmart	1 μ L
DDW	to 10 μ L

Table 5

Digested pRBI2s-single-T7-STS is verified by electrophoresis which is showed below.




Figure 2

It is inferred from the gel electrophoresis that the length of the sample is as expected. The construction is preliminarily correct. We sequenced the constructed plasmid. The result shows that the sequence is the same as expected, indicating that our target plasmid has been constructed successfully.