

CONNECTION

ABSTRACT

This protocol is used to connect two or more pieces of DNA together

BEFORE STARTING

Set up a small box with ice, put DNA and enzymes on it.

Prepare the water bath to 37°C to have Gibson assembly.

1. Select the appropriate connection method according to the experimental situation.

1.1Gibson connection

1.1.1 C112 Gibson connection

This protocol is based on C112 ClonExpress® II One Step Cloning Kit by Vazyme. This protocol is used to connect two pieces of DNA

Preparation of linearized vectors

2. Select an appropriate cloning site on the vector that will be linearized.

3. Vector linearization: the linearized vector can be obtained by digesting the circular vector with restriction enzymes or by reverse PCR.

PCR of the inserts DNA fragments

4. Amplify the insert DNA fragments with homologous sequences (for homologous recombination) of vector-upstream or -downstream by PCR using high fidelity DNA polymerase.

Calculate amount and ratio of linearized vectors and Inserts

5. Detect DNA concentration of linearized vectors and inserts by Nanodrop.

6. Calculation of the amount of vectors:

The optimal amount of vector for the recombination with ClonExpress II is 0.03 pmol, while the optimal amount of insert is 0.06 pmol (molar ratio of vector to insertion is 1:2), as roughly calculated as follows:

The optimal mass of vector = $[0.02 \times \text{number of base pairs}] \text{ ng}$ (0.03 pmol)

The optimal mass of insert = $[0.04 \times \text{number of base pairs}] \text{ ng}$ (0.06 pmol)

Recombination

7. Dilute linearized vectors and inserts before recombination to make sure the loading accuracy. The volume of each component loaded should be no less than 1 μl .

8. Set up the following reaction on ice:

	A	B
1	Linearized Vectors	X μl
2	Inserts	Y μl
3	5 \times CE II Buffer	4 μl

4	Exnase II	2 μ l
5	ddH ₂ O	Add to 20 μ l

9. Use the palm centrifuge to mix the solution in PCR tube.

10. Incubate at 37°C for 30 min and immediately place the tube at 4°C or on ice.

1.1.2 C115 Gibson connection

This protocol is based on C115 ClonExpress® II One Step Cloning Kit by Vazyme. This protocol is used to connect three or more pieces of DNA

Preparation of linearized vectors

2. Select an appropriate cloning site on the vector that will be linearized.

3. Vector linearization: the linearized vector can be obtained by digesting the circular vector with restriction enzymes or by reverse PCR.

PCR of the inserts DNA fragments

4. Amplify the insert DNA fragments with homologous sequences (for homologous recombination) of vector-upstream or -downstream by PCR using high fidelity DNA polymerase.

Calculate amount and ratio of linearized vectors and Inserts

5. Detect DNA concentration of linearized vectors and inserts by Nanodrop.

6. Calculation of the amount of vectors:

The optimal amount of vector for the recombination with ClonExpress II is 0.03 pmol, and the optimal amount of insert is 0.03 pmol (molar ratio of vector to insertion is 1:1), as roughly calculated as follows:

The optimal mass of vector = [0.02 \times number of base pairs] ng (0.03 pmol)

The optimal mass of insert = [0.02 \times number of base pairs] ng (0.03 pmol)

Recombination

7. Dilute linearized vectors and inserts before recombination to make sure the loading accuracy. The volume of each component loaded should be no less than 1 μ l.

8. Set up the following reaction on ice:

	A	B
1	Linearized Vectors	X μ l
2	Inserts	Y ₁ +Y ₂Y _n μ l
3	2 \times clonExpress Mix	4 μ l
4	ddH ₂ O	Add to 20 μ l

9. Use the palm centrifuge to mix the solution in PCR tube.

10. Incubate at 55°C for 15 min and immediately place the tube at 4°C or on ice.

1.2 Goldengate connection

Preparation of linearized vectors

2. Select an appropriate cloning site on the vector that will be linearized.
3. Vector linearization: the linearized vector can be obtained by digesting the circular vector with restriction enzymes or by reverse PCR.

PCR of the inserts DNA fragments

4. Amplify the insert DNA fragments with homologous sequences of vector-upstream or -downstream by PCR using high fidelity DNA polymerase.

Recombination

5. Dilute linearized vectors and inserts before recombination to make sure the loading accuracy. The volume of each component loaded should be no less than 1 μ l.
6. Set up the following reaction on ice:

	A	B
1	Linearized Vectors	10ng
2	Inserts	Y1+Y2.....Yn μ l
3	BSA Enzyme	0.5 μ l
4	BSA Buffer(1g/L)	1 μ l
5	T4 DNA ligase	0.2 μ l
6	10 \times Ligase Buffer	1 μ l
7	ddH ₂ O	Add to 10 μ l

9. Use the palm centrifuge to mix the solution in PCR tube.
10. Program the thermocycler as follows:

	A	B
1	Temperature	Time
2	37°C	30 min
3	37°C	3 min
4	16°C	3 min
5	16°C	20 min
6	37°C	20 min
7	55°C	15 min
8	80°C	15 min

immediately place the tube at 4°C or on ice.

1.3 T4 connection

Preparation of linearized vectors

2. Select an appropriate cloning site on the vector that will be linearized.
3. Vector linearization: the linearized vector can be obtained by digesting the circular vector with restriction enzymes

Enzymatic digestion of DNA fragments

4. Select appropriate enzyme for digestion

Recombination

5. Dilute linearized vectors and inserts before recombination to make sure the loading accuracy. The volume of each component loaded should be no less than 1 μl .

6. Set up the following reaction on ice:

	A	B
1	Linearized Vectors	0.03pmol
2	Inserts	0.3pmol
3	T4 DNA Ligase	0.2 μl
4	10 \times Ligase Buffer	1 μl
5	ddH ₂ O	Add to 10 μl

9. Use the palm centrifuge to mix the solution in PCR tube.

10. Incubate at 16°C for 30 min and immediately place the tube at 4°C or on ice.