

## CONNECTION

### ABSTRACT

This protocol is used to connect two or more pieces of DNA together

### BEFORE STARTING

- Set up a small box with ice, put DNA and enzymes on it.
- Prepare the water bath or mental bath.

Select the appropriate connection method according to the experimental situation.

#### 1.1 Gibson connection

##### 1.1.1 C112 Gibson connection

This protocol is based on C112 ClonExpress® II One Step Cloning Kit by Vazyme. This protocol is used to connect two pieces of DNA.

##### Preparation of linearized vectors

2. Select an appropriate cloning site on the vector that will be linearized.

3. Vector linearization: the linearized vector can be obtained by digesting the circular vector with restriction enzymes or by reverse PCR.

##### PCR of the inserts DNA fragments

4. Amplify the insert DNA fragments with homologous sequences (for homologous recombination) of vector-upstream or -downstream by PCR using high fidelity DNA polymerase.

##### Calculate amount and ratio of linearized vectors and inserts

5. Detect DNA concentration of linearized vectors and inserts by Nanodrop.

6. Calculation of the amount of vectors:

The optimal amount of vector for the recombination with ClonExpress II is 0.03 pmol, while the optimal amount of insert is 0.06 pmol (molar ratio of vector to insertion is 1:2), as roughly calculated as follows:

The optimal mass of vector =  $[0.02 \times \text{number of base pairs}] \text{ ng}$  (0.03 pmol)

The optimal mass of insert =  $[0.04 \times \text{number of base pairs}] \text{ ng}$  (0.06 pmol)

##### Recombination

7. Dilute linearized vectors and inserts before recombination to make sure the loading accuracy. The volume of each component loaded should be no less than 1  $\mu\text{l}$ .

8. Set up the following reaction on ice:

	A	B
1	Linearized Vectors	X $\mu\text{l}$
2	Inserts	Y $\mu\text{l}$
3	5 $\times$ CE II Buffer	4 $\mu\text{l}$

4	Exnase II	2 $\mu$ l
5	ddH <sub>2</sub> O	Add to 20 $\mu$ l

9. Use the palm centrifuge to mix the solution in PCR tube.

10. Incubate at 37°C for 30 min and immediately place the tube at 4°C or on ice.

### 1.1.2 C115 Gibson connection

This protocol is based on C115 ClonExpress® II One Step Cloning Kit by Vazyme. This protocol is used to connect three or more pieces of DNA

#### Preparation of linearized vectors

2. Select an appropriate cloning site on the vector that will be linearized.

3. Vector linearization: the linearized vector can be obtained by digesting the circular vector with restriction enzymes or by reverse PCR.

#### PCR of the inserts DNA fragments

4. Amplify the insert DNA fragments with homologous sequences (for homologous recombination) of vector-upstream or -downstream by PCR using high fidelity DNA polymerase.

#### Calculate amount and ratio of linearized vectors and inserts

5. Detect DNA concentration of linearized vectors and inserts by Nanodrop.

6. Calculation of the amount of vectors:

The optimal amount of vector for the recombination with ClonExpress II is 0.03 pmol, and the optimal amount of insert is 0.03 pmol (molar ratio of vector to insertion is 1:1), as roughly calculated as follows:

The optimal mass of vector = [0.02 × number of base pairs] ng (0.03 pmol)

The optimal mass of insert = [0.02 × number of base pairs] ng (0.03 pmol)

#### Recombination

7. Dilute linearized vectors and inserts before recombination to make sure the loading accuracy. The volume of each component loaded should be no less than 1  $\mu$ l.

8. Set up the following reaction on ice:

	A	B
1	Linearized Vectors	X $\mu$ l
2	Inserts	Y <sub>1</sub> +Y <sub>2</sub> .....Y <sub>n</sub> $\mu$ l
3	2 × clonExpress Mix	4 $\mu$ l
4	ddH <sub>2</sub> O	Add to 20 $\mu$ l

9. Use the palm centrifuge to mix the solution in PCR tube.

10. Incubate at 55°C for 15 min and immediately place the tube at 4°C or on ice.

### 1.2 Golden gate connection

#### Preparation of linearized vectors

2. Select an appropriate cloning site on the vector that will be linearized.
3. Vector linearization: the linearized vector can be obtained by digesting the circular vector with restriction enzymes or by reverse PCR.

#### PCR of the inserts DNA fragments

4. Amplify the insert DNA fragments with homologous sequences of vector-upstream or -downstream by PCR using high fidelity DNA polymerase.

#### Recombination

5. Dilute linearized vectors and inserts before recombination to make sure the loading accuracy. The volume of each component loaded should be no less than 1  $\mu$ l.
6. Set up the following reaction on ice:

	A	B
1	Linearized Vectors	10 ng
2	Inserts	Y1+Y2..... Yn $\mu$ l
3	BSA Enzyme	0.5 $\mu$ l
4	BSA Buffer(1 g/L)	1 $\mu$ l
5	T4 DNA ligase	0.2 $\mu$ l
6	10 $\times$ Ligase Buffer	1 $\mu$ l
7	ddH <sub>2</sub> O	Add to 10 $\mu$ l

9. Use the palm centrifuge to mix the solution in PCR tube.
10. Program the thermocycler as follows:

	A	B
1	Temperature	Time
2	37°C	30 min
3	37°C	3 min
4	16°C	3 min
5	16°C	20 min
6	37°C	20 min
7	55°C	15 min
8	80°C	15 min

immediately place the tube at 4°C or on ice.

### 1.3 T4 connection

#### Preparation of linearized vectors

2. Select an appropriate cloning site on the vector that will be linearized.
3. Vector linearization: the linearized vector can be obtained by digesting the circular vector with restriction enzymes.

#### Enzymatic digestion of DNA fragments

4. Select appropriate enzyme for digestion

**Recombination**

5. Dilute linearized vectors and inserts before recombination to make sure the loading accuracy. The volume of each component loaded should be no less than 1  $\mu$ l.

6. Set up the following reaction on ice:

	A	B
1	Linearized Vectors	0.03 pmol
2	Inserts	0.3 pmol
3	T4 DNA Ligase	0.2 $\mu$ l
4	10 $\times$ Ligase Buffer	1 $\mu$ l
5	ddH <sub>2</sub> O	Add to 10 $\mu$ l

9. Use the palm centrifuge to mix the solution in PCR tube.

10. Incubate at 16°C for 30 min and immediately place the tube at 4°C or on ice.