

1. Obtain *hmas*S mutant fragment by PCR

Adjust the concentration of Mn^{2+} and Mg^{2+} to determine the appropriate mutation rate

PCR system (50 μ L)		PCR system (50 μ L)	
Scpa1-F	2 μ L	Scpa1-F	2 μ L
Hmas-R	2 μ L	Hmas-R	2 μ L
gy9s(scpa1)-177-hmas3	2 ng	gy9s(scpa1)-177-hmas3	2 ng
2 x Mix	25 μ L	2 x Mix	25 μ L
Mn^{2+} (1 mM)	2 μ L	Mn^{2+} (1 mM)	0.5 μ L
Mg^{2+} (25 mM)	2.5 μ L	Mg^{2+} (25 mM)	2.5 μ L
dNTP	1 μ L	dNTP	1 μ L
DDW	14.5 μ L	DDW	16 μ L

PCR system (50 μ L)		PCR system (50 μ L)	
Scpa1-F	2 μ L	Scpa1-F	2 μ L
Hmas-R	2 μ L	Hmas-R	2 μ L
gy9s(scpa1)-177-hmas3	2 ng	gy9s(scpa1)-177-hmas3	2 ng
2 x Mix	25 μ L	2xMix	25 μ L
Mn^{2+} (1 mM)	1 μ L	Mn^{2+} (1 mM)	1.5 μ L
Mg^{2+} (25 mM)	2.5 μ L	Mg^{2+} (25 mM)	2.5 μ L
dNTP	1 μ L	dNTP	1 μ L
DDW	15.5 μ L	DDW	15 μ L

The target fragment is located on the Pccd-K-scpa1 plasmid, and the amplified target fragment is obtained by PCR.

The PCR products are detected by agarose gel electrophoresis, and the correct target fragment is 1248bp. We obtained the correct target fragment, and the sample is purified.

1. Obtain Bio177 vector fragment by PCR

The target fragment is located on the pYB1a-Bio177-Str plasmid, and the amplified target fragment is obtained by PCR. The PCR system and procedure are showed below.

PCR system (50 μ L)	
Bio177-F	2 μ L
Bio177-R	2 μ L
gy9s(scpa1)-177-hmaS3	2 ng
2 x Mix	25 μ L
DDW	20 μ L

The PCR products are detected by agarose gel electrophoresis, and the correct target fragment is 5089bp. We obtained the correct target fragment, and the sample is purified.

2. Golden gate assembly

The *hmaS* mutant fragments and Bio177 vector fragments are connected by Golden gate assembly method, and the system is as follows.

Golden gate system (50 μ L)		Golden gate	
hmaS mut	1 μ L	37 $^{\circ}$ C	60 min
Bio177	1 μ L	55 $^{\circ}$ C	15 min
T4 Buffer	1 μ L	80 $^{\circ}$ C	15 min
rCutsmart	1 μ L		
T4 ligase	0.2 μ L		
BsaI	0.5 μ L		
DDW	5.3 μ L		


3. Colony PCR

After cultured at 37 $^{\circ}$ C for 12 hours, 30 colonies were selected on the plate. The colony PCR system and procedure were as follows.

PCR system (50 μ L)	
Bio177-F	2 μ L


Bio177-R	2 μ L
gy9s(scpa1)-177-hmas3	2 ng
2 x Mix	25 μ L
DDW	20 μ L

The PCR products were detected by agarose gel electrophoresis, and the results were as follows.



4. Sequencing for ensured appropriate mutation rate


Group 1: The final concentration of Mn^{2+} is 0.04 mM




Group 2: The final concentration of Mn^{2+} is 0.03 mM



Group 3: The final concentration of Mn^{2+} is 0.02 mM





Group 4: The final concentration of Mn^{2+} is 0.01 mM

