

8.19

ARAC total 20 μ l, 1.5ml Centrifuge tube

① ddH₂O: 6 μ l

② Buffer: 2 μ l

③ substrate: 10 μ l

KpnI, SalI each 1 μ l

Flick to mix, instantaneous centrifuge

Preparation for sample addition 10min

TETR total 20 μ l

① ddH₂O: 6 μ l

② Buffer: 2 μ l

③ substrate: 10 μ l

PstI, SpeI each 1 μ l

37°C 15 min

80°C 20 min

Glue preparation: 0.5g + 5ml Tae + melt \rightarrow 5 μ l Dye

5 μ l Two kinds of: 5000, 15000 + substrate 20 μ l + loading buffer 4 μ l (10min)

Power: voltage max + time 26min

8.21

Results : tetr with two bands (more than 3000, more than 700 arac with one band 2000-3000)

Note: First turn on the gel imager next to the computer, and then start the software. (When the gel imager is turned on, there will be two yellow lights, "Power" and "Slow" lights.)

Cut glue: open the bottom cabinet of the imager, lay plastic wrap, insert the UV baffle, place the glue in the center, open the TRANS UV knife to cut the required strips, remove the empty glue as far as possible, chop it, transfer it into the empty centrifugal tube, turn off the UV and power supply, remove the plastic wrap and baffle, and wipe the table

The amount of tetr is too small and may not be recycled

tetr plasmid was placed in two 50 μ l tubes at -20°C

arac 0.31g tetr 0.57g

PN 310 μ l PN 570 μ l

PW plus anhydrous ethanol 60ml

(If the volume is too large, add a collection tube) (speed 12000rpm)

1. Column equilibrium, adsorption column 500 μ IBL, 1min
2. Sol, 310 μ l/570 μ l, 50°C
3. Transfer to the adsorption column, place at room temperature for 2min, then transfer for 1min
4. Add 600 μ IPW for 1min
5. Repeat 4
6. Idle for 1 minute, then remove the lid and leave to dry for a few minutes
7. >30 μ IEB (be sure to drop on the white cylinder), replace the tube, and rotate for 2 minutes

After cutting and recovering, the concentration was: arac: 18ng/ μ l, tetr: 15.2ng/ μ l
14.4ng/ μ l

Mixed with 0.50g+50mlTAE+5 μ l dye

digestion

Single enzyme digestion:

A1 arac 18ng/ μ l

arac 44.4 μ l——800ng

Buffer 1.6 μ l (2*0.8)

Enzyme 0.8 μ l (1*0.8)

Double digestion:

A2 Arac (250ng/ μ l) 4 μ l——1000ng

Buffer 2 μ l

Enzyme KpnI+SaII 2 μ l

H2O 12 μ l

Detection

A1:5.88 μ l 5.9 μ l

A2: 2 μ l

Run the glue (120V 400mA 30min) :

Electrophoresis1 Marker

Electrophoresis2 plasmids

Electrophoresis3 A1

Electrophoresis4 A2

8.22

Digestion tetr (T) 107ng/ μ l

Run the glue (120V 400mA 24min) :

Electrophoresis1 Marker

Electrophoresis2 T

Electrophoresis3 A2(finalizing Single enzyme digestion)

Gum cutting: tetr 0.08g arac0.09g

PN 80 μ l

90 μ l

Metallic bath 2min

Rubber recovery

Measuring concentration:

tetr recovery concentration: 5.5ng/ μ l

Arac: 10.9ng/ μ l (A full ng takes 91.74 μ l)

Run the glue A1(120V 400mA 30min)

Loading buffer:Each hole has about 2 μ l

A1:Each hole has about 8 μ l

Electrophoresis1 2 μ l loading buffer

Electrophoresis2 8 μ l A1

Electrophoresis3 marker GL5000

Second electrophoresis:

Electrophoresis1 4 μ l loading buffer 10 μ lA1

Electrophoresis2 2 μ l loadingbuffer 8 μ l A1

Electrophoresis3 marker GL5000

8.23

A1 Cut glue (with about 30ml TAE liquid) and recycle glue

Empty pipe: 1.15g

	Pipe1	pipe2	pipe3	pipe①	pipe②	pipe③
	0.72g	0.57g	0.54g	0.37g	0.64g	0.57g
PN	720	570	540	370	640	570

Glue recovery, eluent EB 50 μ l, centrifuge 3 times

Concentration 7.7ng/ μ l , $n \cdot M = c \cdot v$,To get to 0.1pmol you need about 84 μ l of solution [1057.1500]

Ptcr99a Double enzyme cutting, running glue

3.2g 100ml LB medium

8.24

Solid to liquid medium

The large vial of liquid culture has been spiked with ampicillin

One colony was selected from ① and ② and added to 7ml medium

8.25

Extract plasmid (3ml bacterial solution plus 70 μ lEB)

concentration

Digestion:

DNA 45 μ l

Buffer 2 μ l

Kpn1,Sa2 2 μ l

Glue 30ml

Run the glue

Electrophoresis1 Cutting plasmid
Electrophoresis2
Electrophoresis3 Cutting plasmid 4 μ l
Electrophoresis4 plasmids 4 μ l
Electrophoresis5 marker 1.5 μ l

Weighing after cutting glue: 0.26g
260 μ l PN
ptrc99a, enzyme digestion 9.4ng/ μ l, need 0.01pmol=0.01*10⁽⁻¹²⁾mol=10⁽⁻¹⁴⁾mol
4159*650*10⁽⁻¹⁴⁾=2.7*10⁽⁻⁸⁾g=27ng
Volume: 27 \div 9.4=2.9 μ l

8.26

The receptive + connected plasmid (arac) was added with 25ngDNA per 50 μ l cell
Yesterday: ptrc99a:27ng, arac:165ng
A total of 192ng(but the arac is increased, it should be greater than 192ng) is calculated at 200ng
For 100 μ l receptive cells (50 μ l cells :25ngDNA), 50ngDNA should be added
200 ng / 2.9 + 84 + 1 + 1 \approx 90 μ l
Take :22.5ml and attach the plasmid
Sterilization: solid medium, liquid medium (solid medium,4.8g powder +150ml pure water, add 0.015gAmp after sterilization
Liquid medium,2.5g powder +100ml pure water)

8.26 PM

100 μ l K12+22.5 μ l connected plasmid
Follow the K12 instructions
arabinose(inducer) was not added before the bacteria were attached to the plate, and arabinose was added after the bacteria were attached
arabinose is powder, 0.03g(150ml medium)
0.03g arabinose+1100ul pure water
Eleven plates were poured and 100 μ l arabinose solution was added to each plate
The green colonies were cultured and screened for successful connection

8.28 PM.

Enzyme digestion of tetr fragments:
t1 16.1ng/ μ l, t2 19.2ng/ μ l
Blp1 2 μ l Buffer

8.30 PM.

Run tetr cut glue cut waste, no strip
Start from the large fragment of tetr again
Arabinose dosage :1mM concentration, 100 μ l per plate
Maracbinose = 342.3 g/mol
Total 5 plates,500 UL solution,5*10⁽⁻⁴⁾mol aracbinose

Mass : $5 \times 10^{-4} \times 342.3 = 0.17\text{g}$

That is, 0.17g + 500µl pure water

The bacteria cultured in 8.26 (placed in the refrigerator) were observed by fluorescence microscope, and the phenomenon was not obvious, and the colonies were in large pieces, which may have been caused by the lack of smooth coating during inoculation. Today, the colonies with slight bright spots in fluorescence observation were selected → 200µl LB + Amp liquid medium dilution → coated to LB plate (the incubator was not placed because it was closed).

In the next experiment, green colonies were found, and each observed one was marked at the bottom of the plate during fluorescence observation. The labeled colonies were picked into liquid medium, and the plasmid could be extracted after culture

8.31

1. Column balance (1.08)
2. Weighing of rubber block: 0.05g (0.1g → 100µl) → Adding 50µl PN, it is difficult to dissolve, and adding 30µl PN (weighing centrifugal tube problem) → resulting in a volume of 160µl
3. Add CA2, 2min, then centrifuge
4. Add 600µl PW
5. Add EB Eluent about 35µl

(Repeat centrifugation twice during the collection phase)

Change 1.5ml centrifuge tube for collection (concentration to be measured)

Attempted digestion 6.9ng

Tetr 19µl	DNA 1µg	13µl
Buffer 2µl	Pst1 1µl	Spe1 1µl
B1P1 1µl	Buffer 2µl	

37°C 1h

The tetr gene has no band and cannot be recycled

~700bp glue recovery

PM: The tetr fragment was recovered and the concentration was measured at 8.6ng/µl

9.11

Enzyme digestion arac, tetr

Electrophoresis:

Electrophoresis1 maker

Electrophoresis2 arac

Electrophoresis3 tetr

Electrophoresis4 tetr

When the plasmid is used up, the plasmid is extracted from the bacterial solution, but the storage time is too long, the bacteria sink to the bottom, not the suspension of the fresh bacterial solution, which may affect the result

Plasmid concentration: a1a2 is the same, t1t2 is the same

A1 493.8ng/μl

A2 521.9ng/μl

T1 149.7ng/μl

T2 119.1ng/μl

9.12

Digestion a1 t1

A1 substrate 2.02μl t1 substrate 6.68μl

Buffer 2μl buffer 2μl

Kpn1 Sa2 2μl Pst1 Spe1 2μl

DdH2O 14μl DdH2O 9.32μl

Electrophoresis:

Electrophoresis1 maker

Electrophoresis2 arac

Electrophoresis3 tetr

Rubber recovery:

Arac 0.16g tetr0.2g

160μl 200μl

8.8ng/μl 4ng/μl

8.8ng/μl 0.5μl

DNA 6μl

Buffer 2μl

BLP1 1μl

DdH2O 11μl

37°C 1.5h

Electrophoresis (30ml) 30min

Two target bands can be seen, but the amount is small and stored at 4°C

arac single enzyme cut, cut glue, stored in the upper layer of 4°C refrigerator

Glue recovery, concentration measurement, 5.8ng/μl

9.14

Pttrc99a 53.3ng/μl

Arac Kpn1 Sa2 493.8ng/μl

Tetr Pst1 Spe1 149.7ng/μl

Pttrac99a Kpn1 Sal1 53.3ng/μl 18.8μl

Maker tetr arac ptrc

① ② ③ ④

Empty pipe: 1.00g

Pttrc99a: 0.16g 16.6 3μl

Tetr: 0.10g

Arac: 0.09g 15.1 12μl

Glue recovery, concentration not measured

9.17

Preparation: Gun tip, PCR tube, sterilization of medium, inverted plate

Connections:

arac 12 μ l

ptrc99a 3 μ l

Buffer 2 μ l

Ligase 1 μ l

H₂O 2 μ l

16:00~18:00, 25°C for 2h, 4°C storage

Tomorrow:

Open the ice machine, prepare the liquid medium, kanamycin, arabinose

Feelings of ice melting

Inducer is added before coating

The connecting system is refrigerator 4°C at the far gate of 307

9.20

50 μ l \rightarrow 25ngDNA

906.4ng/ μ l 100 μ l \rightarrow 50ngDNA

Dilution :0.5 μ l of the connection product is added to 4.5 μ l ddH₂O

1 μ l connection product \rightarrow 100 μ l IK12

Coating plate 5:

IPTG+ Arabinose :2 plates

IPTG:1 board

Arabinose :1 plate

None :1 board

(1mM arabinose :3.4g \rightarrow 10ml water

1mM IPTG 0.002383g/mol)

0.1M arabinose :0.15g \rightarrow 10ml water

0.1M IPTG:0.238g \rightarrow 10ml water

Memo No. ()

Date

ARAC

(第: 20 ml)

1.5 ml

① ddH₂O 6 ml

② Buffer 2 ml

③ 抗体 10 μl

(KpnI, SalI 各 1 μl) 加样: 10 min

轻轻混匀: 轻轻混匀

20 ml TE TR

① ddH₂O 6 ml

② Buffer 2 μl

③ 抗体 10 μl

Post I - SpelI 各 1 μl 60°C 20 min

凝胶: 0.5 g + 50 ml TAE + 酶消化 → 5 ml 试剂

5 ml 两种 5000, 15000

★ 抗体 20 ml + loading buffer

功: 电泳: 10 min

时: 10 min

V. 10 min

(准备1支) (准备2支)
准备: 2支吸附枪, 4支离心管

Memo No. _____
Date / /

胶回收

2支吸附枪收集管 ^{高加} 500 μ l BL \rightarrow 12000 / 1min 倒废液
切胶 \rightarrow 离心管 \rightarrow 称量 (araC, tetr 各一个)
加PN \downarrow 0.1g 胶 / 100 μ l
金属浴控温水浴锅 50°C 至全溶 (大于 800 μ l 取各针吸附枪分批加)
 \rightarrow 2支吸附枪, 各加一个胶溶液 \leftarrow \rightarrow 放 2min (空) \rightarrow 12000 / 60s

PW 600 μ l + 无水乙醇 . 12000 / 60s . 倒废液

\downarrow
12000 / 2min . 室温放置 5min . 开盖

\downarrow
换新离心管 EB ? μ l 滴在枪头
 \downarrow 擦胶追加. 问老师

室温放置 2min . 12000 / 2min . 离心管中为收集的DNA溶液

制冰: 厚担开. 水笼又供水

金属浴角加热: 调 time, temp \rightarrow start 手机计时

(仪器也不准)

8.21 下午 1. 5000 marker 2. tetr arac 3. Memo No. _____ Date / /

酶切电泳图拍照 结果: tetr 两条带 (3000多, 700多) arac 一条带 2000~3000

√ 注意: 先开电脑旁的凝胶成像仪, 再启动软件, (凝胶成像仪
开启会有两个黄个灯亮, Power 和 Slow 亮)

- 切胶 开成像仪底柜, 铺保鲜膜 插紫外挡板
- 胶放置正中 开 TRANS UV 刀切所需条带
- 尽量剔除空胶 切碎 转移入空离心管
- 关 UV、电源
- 去保鲜膜、挡板
- 擦净台面

tetr 量过少 可能无法回收

提 tetr 质粒 两管 50μl 置 307 -20℃

arac - 0.51g | tetr 0.57g
 PN 310 ml | PN 570 ml

Δ PW 加无水乙醇 60ml

1. 柱平衡 (多加一个吸附柱 → 收集管) 12000rpm

2. 溶解 310 μl / 570 μl 50°C 1min

3. 分别转入吸附柱室温放 2min 后 1min

4. 600 μl PW 1min

5. 重复 1min

6. 空转, 后开盖晾几分钟 2min

7. >30 μl EB (一定要滴到白色柱体上) 换新管 2min

切割回收后浓度: arac: 18 ng/μl | tetr: 15.2 ng/μl

切割 0.5g + 50 ml TAE + 5 μl 染料

• 切割 arac 18 ng/μl 稀释 5 倍

arac 44.4 μl — 800 ng

buffer 1.6 ml (2 × 0.8)

酶 XhoI 0.8 ml (1 × 0.8)

arac (250 ng/μl) 4 μl — 100 ng

buffer 2 μl

酶 KpnI + SalI 2 μl

水 12 μl

稀释 100

20 μl

检测 (跑胶)

a提: 5.88 ml 5.9 ml

a取: 2 ml

泳 - = 三

Marker 泳 - a取 a取

120V 400mA 30min

8.22 下午 107 ng/ml

酶切 tetr 提 最后 (要做单酶切)

泳 - = 三 四 (的)

Marker 提 a取

24min

切胶 tetr 0.08g arac 0.09g

DN 80 ml 90 ml 金属浴 2 min

切胶回收

测浓度:

tetr 回收浓度: 5.5 ng/ml

arac: 10.9 ng/ml (满 1ng 要 91.74 ml)

8.22晚

① 对α提跑胶. ~~切胶. 胶回中~~

loading buffer: 每个孔大约有2μl.

α提: 每个孔大约有8μl.

120V. 400mA. 30min.

泳道一

2μl loading buffer
8μl α提

Marker.
GL5000

第二次电泳

GL5000

泳道一

每个孔有4μl loading buffer.
10μl 左右 α提.

2μl β-buffer
8μl α提

Marker. 2.2

2.2

Marker

2.2

2.2

2.2

2.2

2.2

2.2

2.2

2.2

2.2

2.2

Memo No.
Date

1.92 1.15 Memo No. _____
Date / /

8.23 下午. → 配 30ml 左右 TAE 液

a 提 切胶 胶回收

空管	1.15g	管 1	管 2	管 3	管 ①	管 ②	管 ③
	0.72g	0.57g	0.54g	0.78g	0.61g	0.57g	0.57g
	0.81g	0.39g	0.83g	0.83g	0.83g	0.83g	0.83g
PN	720	570	540	780	610	570	570

胶回收, 洗脱液 EB 50ml. 离心 3次

浓度 7.7 ng/ml. $n \cdot M = c \cdot V$ 到 0.1 pmol 需 84ul 左右溶液

[1057.1500]

ptrc99a 双酶切: 跑胶 = = =

8.23 晚上

3.20g 100ml LB 培养基

8.24 下午

固 → 液

大瓶液体培养基 经 10 倍稀释

从 1, 2 分别挑取一个菌落 加入到 7ml 培养基

8.26 上午

感受态 + 连好的质粒 (arac) 每 50 μ l 细胞加 25 ng DNA

昨天算 pUC 99A: 27 ng arac: 165 ng

共 192 ng (但 arac 加多了, 应大于 192 ng) 下按 200 ng 计算

100 μ l 感受态细胞 (50 μ l 细胞: 25 ng DNA), 应加 50 ng DNA

\rightarrow 200 ng / (29 + 84 + 1 + 1) \approx 90 μ l

取 22.5 μ l 连好质粒

灭菌: 固培: 液培 (固体培养基 4.8g 粉末 + 150 ml 纯水, 灭菌后加 0.015g Amp)

液体培养基 2.5g 粉末 + 100 ml 纯水

8.26 下午

100 μ l K12 + 22.5 μ l 连好质粒

按说明书做

★ 菌接到平板前未加 arabinose (诱导物), 之后补加的 arabinose

arabinose 是粉末 (取 0.03g (培养基 150 ml))

0.03g arabinose + 1100 μ l 纯水

倒了十一个平板, 每个平板加 100 μ l arabinose 溶液

培养, 筛长出绿色菌落为连接成功

8.28 下午

酶切 tet^r 片段: ~~PstI~~ ~~BlnI~~

t₁ 16.1 ng/ μ l t₂ 19.2 ng/ μ l

BlnI. 2 μ l Buffer.

酶 1 μ l.

8.30 晚.

跑 tetr 彻底，无杂带 (control 铁质培养基 + 奈曼液)
重新从 tetr 大片段开始。 p₁₁ 5 = 2pp 20g 第一天

阿拉伯糖加量 (1 mM 浓度: 每平板 100 μl 阿拉伯糖 14.00)

Maracbinose = 342.3 g/mol

共 5 个板 500 μl 溶液 500 · 5 × 10⁻⁴ mol arabinose

质量 = 5 × 10⁻⁴ × 342.3 = 0.171 g (约) 新液: 新固: 菌液

即 0.17 g + 500 μl 纯水 (1 mM) + 奈曼 p₂₁₅ 培养基平板

荧光显微镜观察 8.26 培养的菌 (放冰箱); 现象不明显
菌落成大片可能是接种时未涂匀, 今天将荧光观察稍有亮
点菌落挑取 → (200 μl) LB + Amp 液体培养基稀释 →
涂至 LB 平板 (未放培养箱因关门了) 未测是 *serenid*

下次实验, 挑绿色菌落, 荧光观察时, 有观察到个, 在平板
底部标记, 将有标记的菌落挑至液体培养基,
培养后可提质粒

8.31

Memo No. / /
Date / /

tert 781 bp 收回

1. 称取 ~~0.24~~ ^{0.24} ~~1.10~~ ^{1.08} g DNA

2. 胶块称重: 0.05g (0.1g → 100μl 1P → 加入 50μl PN
溶液和 20μl PN (称量管问题)

⇒ 等效体积 160μl

3. 加入 CA2. 室温 2min

4. 加入 600μl PN

5. 加入 EB 洗液 约 35μl

(收集到较重条带 2次)

→ 换 1.5ml 离心管收集 (浓度检测)

6.9g → 6.9g

走试酶切 6.9ng

totr 1μl

DNA 1μg 13μl

Buffer 2μl

PstI 1μl SpeI 1μl

BIPI 1μl

Buffer 2μl

37°C 1h

tert 基因无条带

无法回收

~700bp 找回

下午: 胶回收长 tert 片段. 测定浓度 8.6 ng/μl

9.11 上午

酶切 arac, tetr. 体系同第1页 (arac-双酶切)

电泳 泳一 泳二 (1.0) 泳三 泳四

Marker arac tetr tetr

9.11 下午

质粒用完了, 从菌液提质粒, 但保藏时间过长, 菌沉淀, 不是

新鲜菌液的悬浊液, 可能影响结果

质粒浓度: a1 a2 是一个东西 t1 t2 是一个东西

a1 493.8 ng/ul

a2 521.9 ng/ul

t1 149.7 ng/ul

t2 191.1 ng/ul

9.12 上午

酶切 a, t1

a. 质物 2.02 ml

buffer 2 ml

KpnI 50U 2 ml

ddH2O 10 ml

t1. 质物 6.68 ml

buffer 2 ml

PvuII 50U 2 ml

ddH2O 10 ml

电泳 右→左

泳一

泳二

泳三

Marker

arac

tetr

5.68 x 17-x

Memo No.
Date

8.8 x 10⁻⁵ x = Memo No.
Date

9.12 下午

切膜回收

Arac 0.16g

160 μ L

8.8 ng/ μ L

8.8 ng/ μ L

DNA

Butter

BIP1

ddH₂O

37°C 1.5h

电泳 (30mL) 30min

NIP 可见两条目标条带但量较少

9.12 晚

Arac 溶解切切膜 存4°C. 冰上解冻. 立即溶解在下
膜回收. 现浓度. 5.8 ng/ μ L.

x-11 x 882
Memo No. 8.8 x 8.8
Date

Memo No.
Date

9.14 18.8

15 51.9

ptro99a 53.3 ng/ml
Jy000 Jy011

51.81 18.8

arae. kpnI. SalI 8.8 493.8 ng/ml

Tetr. PstI SpeI. 149.7 / ml 8.8

ptro99a. kpnI. SalI. 53.3 ng/ml 18.8 μl

Multer Tetr. Arae 19.7

① ② ③ ④

$\frac{1.00g}{\sqrt{2}} = 0.707g$

↓

ptro99a 0.16g 412

tetr: 0.10g

arae: 0.09g 15.1 12 μl

胶回收 12 μl 18.8

9.20 晚

稀释 100ml 100ml → 50ng DNA

50ml → 25ng DNA

906.4ng/1ul 100ml → 50ng DNA

稀释: 0.5ml 直接产物 加入到 4ml 4.5ml ddH₂O 中

1ul 连接产物 → 100ul K12

涂板 5个:

IPTG + 阿拉伯糖: 1板

IPTG: 1板

阿拉伯糖: 1板

无: 1板

1mM 阿拉伯糖: 3.4g → 10ml 水

1mM IPTG: 0.002383 mol/g/ml · 10⁻² × 10³ mol

238.3g/L = 1M

0.238g/L = 1mM 10ml

0.238g

0.1M 阿拉伯糖: 0.15g → 10ml 水

0.1M IPTG: 0.238g → 10ml 水