

## PCR product purification

### ABSTRACT

This protocol is used to purify PCR products (in case of single band detected by gel electrophoresis)

### BEFORE STARTING

Prepare CP Buffer and Wash Buffer reagents

1. Add 4-5 times CP Buffer to the PCR product, Blow and mix well.
2. Insert a HiBind® DNA Mini Column into a 2 mL Collection Tube.
3. Inhale the mixed liquid into the HiBind® DNA Mini Column
4. Centrifuge at 10000x g speed for 1 minute.  
10000 x g, Room temperature, 00:01:00
4. Discard the filtrate and reuse the collection tube.
5. Add 600 µL Wash Buffer
7. Centrifuge at maximum speed for 1 minute.  
15000 x g, Room temperature, 00:01:00
8. Discard the filtrate and reuse collection tube.
9. Repeat step 6~8 once.
10. Centrifuge the empty HiBind® DNA Mini Column for 2 minutes at maximum speed to dry the column matrix.
11. Transfer the HiBind® DNA Mini Column to a clean 1.5 mL microcentrifuge tube.
12. Add 30-100 µL Elution Buffer or sterile deionized water directly to the center of the column membrane.  
The efficiency of eluting DNA from the HiBind® DNA Mini Column is dependent on pH. If using sterile deionized water, make sure that the pH is around 8.5.
13. Let sit at room temperature for 1 minute.
14. Centrifuge at maximum speed for 1 minute.  
15000 x g, Room temperature, 00:01:00  
This represents approximately 70% of bound DNA. An optional second elution will yield any residual DNA, though at a lower concentration.
15. Suck out the solution from the tube and re-add it to the center of the column membrane to give a second centrifuge.  
15000 x g, Room temperature, 00:01:00
16. Test the concentration and purity of DNA using NanoDrop.
17. Store DNA at -20°C.