

2024 NEFU_China

1. Take pLB1s-ndmDtB as template (template plus 20~50ng/50 μ L), use high-fidelity enzyme, concatenate primers for PCR (25 cycles, reduce unnecessary mutations), amplify 3~4 tubes.
2. Purify the PCR product and do not recover the gel.
3. Take 200~300ng for DpnI digestion.
4. Transfer to BW/TPR sensory state, A+S resistance, incubate at 37 $^{\circ}$ C for 12~16h, try to cultivate until the colony is bigger.
5. Prepare a deep well plate, add 600 μ L of induced transformation medium (LB+0.5 mM theophylline +A100+S50+0.2 % arabinose+1 mM IPTG) into each well in three shots.
6. Pick the monoclones into each well (select large and full colonies as much as possible)
Note: choose four wells to pick the unmutated TPR+DtB (prepared in advance) as the control experiment, the other 92 wells to pick the TPR+DtBs.
7. Cultivate in 96-well plate shaker at 30 $^{\circ}$ C and 850rpm for 24h.
8. Take 200 μ L in ultra-clean and put it into 96-well plate, and detect it under the enzyme labeling instrument.
9. Record the group with stronger fluorescence than the control group, take 50 μ L of bacteria, and 300 μ L of sample for HPLC detection (400 μ L of bacterial solution left in the deep-well plate in the previous experiment was used for this step).
10. Sequence the mutants with improved yield in the HPLC assay to identify their mutation sites. Note: This experiment is applicable to TPR-mCherry, which is a sensor with very slow colony fluorescence formation. If a sensor with fast and strong fluorescence is developed later, the first step is to perform a preliminary screening of the naked eye color to reduce the workload.

Note: This experiment is applicable to TPR-mCherry, which is a sensor with very slow colony fluorescence formation. If a sensor with fast and strong fluorescence is developed later, the first step is to perform a preliminary screening of the naked eye color to reduce the workload.