

Time: 2024.08.21-2024.08.25

1. Experiment: Co-culture of NK92 and human hepatoma organoids

2. Time: 2024.08.21-2024.08.25

3. Member: Xiaoyuan Chen, Hanyue Liu

4. Material:

(1) Organoids: human hepatoma organoids

(2) Reagent: 10 μ M Y-27632, 0.01 μ M Prostaglandin E2 (PGE2; R&D Systems)

5. Method:

For organoid cytotoxicity assays, organoids were seeded on a Matrigel layer and incubated with the respective NK92 cells for 8 h as follows.

- (1) 48-well standard culture plates with a growth area of 1 cm² (Greiner Bio-One) were first moistened using culture medium.
- (2) Subsequently, each well was evenly covered with 35 μ L undiluted Matrigel which was allowed to solidify overnight at RT.
- (3) Confluent organoids were collected, mechanically sheared, pelleted, washed, and seeded at a split ratio of 1:2.5. Organoids were resuspended in 150 μ L of the respective culture medium supplemented with 10 μ M of Y-27632 per assay replicate.
- (4) The organoid suspension was carefully added to the center of the Matrigel-covered wells, respectively. Organoids were grown for 24 h before supplement of NK92 cells in 500 μ L of medium without Y-27632.
- (5) For standard cytotoxicity assays, $\sim 10^5$ organoids cells were seeded per well of the 48-well plate (1 cm²), which was defined as organoid density of 100%. NK92 cells and their CAR-engineered derivatives were pelleted, washed, and resuspended in co-culture medium as indicated below.
- (6) Cells were counted using a hemocytometer, and the required number of cells in a total volume of 500 μ L (or 1 mL for long-term co-cultures) per well of a 48-well plate was co-incubated with target cells for 8 h at 37 °C.
- (7) The following co-culture media were used: For human organoids, the medium contained complete normal or tumor medium lacking nicotinamide. 0.01 μ M Prostaglandin E2 (PGE2; R&D Systems) was added to induce a cystic phenotype.