

Time: 2024.08.17-2024.08.20

1. **Experiment:** Human hepatoma organoids culture and passage
2. **Time:** 2024.08.17-2024.08.20
3. **Member:** Xiaoyuan Chen, Hanyue Liu
4. **Material:**

| Name | Supplier or Formulation |
|-------------------------|--|
| DMEM basic | Invitrogen |
| Stopping culture medium | DMEM basic, Penicillin-streptomycin (100 µg/mL), FBS (10%) |
| Complete culture medium | Advanced DMEM/F12, Penicillin-streptomycin (100 µg/mL), Double-antibody (1×), HEPES (10 mM), GlutaMAX-I (1×), A83-01 (500 nM), Y-27632 (10 µM), N-Acetylcysteine (1.56 mM), Nicotinamide (10 mM), FGF10 (10 ng/mL), B27 supplement (1×), Forskolin (10 µM), Wnt3A conditioned medium (30%), R-spondin conditioned medium (2%), Noggin conditioned medium (4%) |
| Ordinary culture medium | Advanced DMEM/F12, Penicillin-streptomycin (100 µg/mL), Double-antibody (1×), HEPES (10 mM), GlutaMAX-I (1×), A83-01 (500 nM), Y-27632 (10 µM), N-Acetylcysteine (1.56 mM), Nicotinamide (10mM), FGF10 (10 ng/mL), B27 supplement (1×), Forskolin (10 µM), EGF (50 ng/mL), Wnt3A conditioned medium (30%), R-spondin conditioned medium (2%), Noggin conditioned medium (4%) |
| Basal culture medium | Advanced DMEM/F12, Penicillin-streptomycin (100 µg/mL), Double-antibody (1×), HEPES (10 mM), GlutaMAX-I (1×) |
| Penicilin-streptomycin | Invitrogen |
| TrypLE | Invitrogen |
| Collagenase, Type 2 | Invitrogen |
| Y-27632 | Selleck |
| Wash medium | DMEM basic, Penicillin-streptomycin (100 ug/mL) |
| Digestion solution | Collagenase, Type 2 2.5 mg/mL Y-27632 10 uM |
| PBS 10X (1000 mL) | Na ₂ HPO ₄ 1.44 g KH ₂ PO ₄ 0.24 g Add distilled water to 800 mL |
| Matrigel | Corning (356231) |

Time: 2024.08.17-2024.08.20

5. Method:

(1) Culture of human Hepatoma organoids

- ① For tumor tissue samples, approximately 1 cm³ tissue blocks are placed in 10 cm culture dishes and washed twice with clean culture medium.
- ② Scissors and forceps were used to mince the tissue samples and transferred into 15 mL centrifuge tubes containing digestion solution (basal culture medium containing 2.5 mg/mL Type II collagenase and 10 μM Y-27632).
- ③ The tubes are placed on a tube rotator in a 37 °C cell culture incubator for thorough mixing and digestion for approximately one hour. Digestion can be stopped when no large tissue particles are visible in the centrifuge tubes.
- ④ The digestion product is filtered through a 70 μM filter and washed for several times.
- ⑤ The filtrate is centrifuged at 1700 rpm for 5 minutes. The supernatant was discarded.
- ⑥ 3-5 mL stopping culture medium is added. The centrifugation step is repeated at 1700 rpm for 5 minutes.
- ⑦ This process is repeated once more with 2 mL complete culture medium (Hepatoma tissue) or ordinary culture medium (normal liver tissue).
- ⑧ Based on the amount of cell precipitation obtained, 150 μL of the respective culture medium and 100 μL of Matrigel (Store on ice during the experiment) are seeded in suspension 6-well plates, or 30 μL of the respective culture medium and 20 μL of Matrigel are seeded in suspension 24-well plates.
- ⑨ After seeding, the plates are left undisturbed in a cell culture incubator (37 °C, 5% CO₂, and humidity maintained).
- ⑩ After 1 hour incubation, the seeded plates are supplemented with liquid medium (3 mL per well for 6-well plates; 1 mL per well for 24-well plates).

(2) Passage of human Hepatoma organoids

- ① Monitor the organoid culture via microscopy and passage the organoids when density reaches approximately 70% in suspension.
- ② Disperse the organoids using a 1 mL pipette gun into plate wells, then transfer into a 15 mL centrifuge tube.
- ③ Spin at 1500 rpm for 5 min (if the matrix gel-containing precipitate is suspended after centrifugation, the centrifuge can be slowed down)
- ④ Discard the supernatant to 1 mL, and blow the precipitate 80-100 times with a 1 mL pipette gun.
- ⑤ Under microscope observation, the organoid should form a mass of 30-50 μM.
- ⑥ Wash the solid organoids with pre-cooled phosphate-buffered saline (PBS) until no visible matrix gel remains following centrifugation. In cases where solid organoids are resistant to dissociation, the supernatant may be removed post-centrifugation, followed by the addition of 1-2 mL of TrypLE for enzymatic digestion at 37 °C within a cell culture incubator for approximately 8-10 minutes.
- ⑦ Add 5 mL of termination medium and centrifuge at 1000 rpm for 3 min to obtain the cell precipitate.
- ⑧ Resuspend the cells with the appropriate amount of corresponding medium and matrix gel, inoculate the cells on the suspension plate, and control the substitution ratio at 1:2~1:3. The pipette gun should not be stuck to the wall of the centrifuge tube to prevent cross contamination.