

Time: 2025.05.10-2025.05.29

- Experiment:** Flow Cytometric Analysis of NK Cell Activity
- Time:** 2025.05.10-2025.05.29
- Member:** Xudong Tang, Yang Jin, Binxuan Zhang, Kaiqing Zhang, Xuanton Liu
- Materials:** Flow cytometer, CD107a-APC, CD56-PE, CD69-FITC, Ice-cold PBS Staining buffer, Brefeldin A, 4% paraformaldehyde (PFA), 0.1% Triton X-100

**5. Method:****(1) Cell preparation:**

NK cells pre-activated by B34G35R, B51G35R, and B51G9 are collected from previously proliferated cells, with  $0.5-1 \times 10^6$  PBMCs or  $0.1-0.5 \times 10^6$  purified cells per condition.


The cells are resuspended in pre-warmed (37 °C) culture medium at appropriate density.

**(2) Surface staining:**

- CD69 expression was analyzed in all twenty groups demonstrating cytotoxic activity using flow cytometry.
- Post-cytotoxicity assay, cells in each well were washed once with 1 mL of ice-cold PBS and resuspended in 100  $\mu$ L staining buffer.
- Cells were stained with CD69-FITC (BioLegend, 310904), CD56-PE (BioLegend, 362508), CD107a-APC (BioLegend, 328620) to identify NK cells.
- Activated NK cells were defined as CD56<sup>+</sup>CD69<sup>+</sup>, while total NK cells were CD56<sup>+</sup>.
- CD69-FITC (1:100), CD107a-APC and CD56-PE (1:50) were added. Samples were incubated at 4 °C for 30 min protected from light. Control tubes were prepared.
- 200  $\mu$ L ice-cold Flow Staining Buffer was added. The tube was centrifuged at  $300 \times g$  for 5 min. And the supernatant was aspirated completely.
- The cells are resuspended with 100-200  $\mu$ L Fixation Buffer, and incubated at room temperature for 15-30 min.
- 200  $\mu$ L Flow Staining Buffer was added and the tube was centrifuged at  $300 \times g$  for 5 min. The supernatant was aspirated completely.

**(3) Acquisition:**

- The suspension was transferred to flow cytometry sample tubes.
- Analysis Strategy:  
Gate to exclude doublets (using FSC-A vs FSC-H).  
Gate to exclude dead cells (using a viability dye).  
Identify the target cell population (NK cells: CD3<sup>-</sup> CD56<sup>+</sup>).  
Analyze the expression levels of CD69 and CD107a on the surface of NK cells.

**6. Result:**

**Figure.1** CD107a and CD69 expression were analyzed in effector cells after stimulated by IL-2 and its mimics. Data are representative of at least three independent experiments (\*\* $P < 0.01$ , \* $P < 0.05$ ).