Time: 2025.05.10-2025.05.29

- 1. Experiment: Flow Cytometric Analysis of NK Cell PFN & GrB
- **2. Time:** 2025.05.10-2025.05.29
- 3. Member: Xudong Tang, Yang Jin, Binxuan Zhang, Kaiqing Zhang, Xuantong Liu
- **4. Materials:** Flow cytometer, Anti-human Perforin antibody (BioLegend, 308104), Anti-Human Granzyme B Antibody (BioLegend, 515404), Ice-cold PBS Staining buffer, Brefeldin A, 4% paraformaldehyde (PFA), 0.1% Triton X-100

5. Method:

(1) Cell preparation:

- (1) NK cells pre-activated by B34G35R, B51G35R, and B51G9 are collected from previously proliferated cells and Co-cultured with target cells (K562) in logarithmic growth phase with effector cells (pre-stimulated NK cells) at a specific effector-to-target ratio in a 37 °C, 5% CO₂ incubator.
- ② Add Brefeldin A (diluted 1:1000) to the co-culture system 12 hours after initiation. Continue culture for a total of 46 hours.
- 3 After co-culture, collect cells into flow cytometry tubes. Wash once with pre-chilled PBS, centrifuge at 300 × g for 5 minutes, and discard the supernatant.
- 4 Resuspend cells in 100 μL flow cytometry staining buffer for cell counting.

(2) Surface Staining

- ① Dosage example for a 100 μL reaction system: CD3-PerCPCy5.5 (1:100), CD56-PE (1:100).
- ② Add antibodies to the cell tube, gently mix by pipetting, and incubate at 4 °C in the dark for 30 minutes.
- ③ Wash with 2 mL staining buffer, centrifuge at 300 × g for 5 minutes, and discard supernatant. Repeat washing once.

(3) Fixation and Membrane Permeabilization

- ① Resuspend cells in 100 μL fixative, incubate at room temperature in the dark for 15-20 minutes.
- 2 Add 2 mL staining buffer, centrifuge, discard supernatant, and thoroughly remove fixative.
- 3 Resuspend cells in 100 μL lysis buffer and incubate at room temperature in the dark for 15-20 minutes.

(4) Intracellular Staining

- ① Dilute intracellular antibodies in lysis buffer: Perforin-FITC (1:20), Granzyme-BV421 (1:20), centrifuge and discard lysis buffer.
- 2 Add the intracellular antibody cocktail to the cell tube, gently pipette to mix, and incubate at room temperature in the dark for 30-45 minutes.
- (3) Wash with 2 mL lysis buffer, centrifuge at 300 × g for 5 minutes, and discard the supernatant. Repeat the wash once.
- 4 Resuspend cells in 300-500 μL flow cytometry staining buffer and analyze immediately.

(5) Data Analysis:

- 1) The suspension was transferred to flow cytometry sample tubes.
- 2 Analysis Strategy:

Gate to exclude doublets (using FSC-A vs FSC-H).

Gate to exclude dead cells (using a viability dye).

Identify the target cell population (NK cells: CD3⁻ CD56⁺).

Analyze Perforin and Granzyme B expression: Percentage of positive cells (%) and Mean Fluorescence Intensity (MFI).

Time: 2025.05.10-2025.05.29

6. Result:




Figure.1 Killing activity of NK cells and T cells pre-stimulated by IL-2 (control), B34G51, B34G35, B51G9. The PFN and GrB expression levels were assessed by flow cytometry. Data are representative of at least three independent experiments (**P < 0.01, *P < 0.05).