Time: 2025.09.10-2025.09.15

- 1. Experiment: Flow cytometric analysis of NK cell PFN & GrB under acidic condition
- **2. Time:** 2025.09.10-2025.09.15
- 3. Member: Xudong Tang, Yang Jin, Binxuan Zhang, Kaiqing Zhang, Xuantong Liu
- **4. Materials:** Flow cytometer, Anti-human Perforin antibody (BioLegend 308104), anti-Human Granzyme B Antibody (BioLegend 515404), Ice-cold PBS Staining buffer, Brefeldin A, 4% paraformaldehyde (PFA), 0.1% Triton X-100

5. Method:

- (1) Cell preparation:
 - ① Co-culture target cells (K562) in logarithmic growth phase with effector cells (pre-stimulated NK cells) at a specific effector-to-target ratio in a 37 °C, 5% CO₂ incubator, under acidic condition of 15 mM lactic acid (pH 6.4), and neutral pH of 7.4 respectively.
 - ② Add Brefeldin A (diluted 1:1000) to the co-culture system 12 hours after initiation. Continue culture for a total of 46 hours.
 - 3 After co-culture, collect cells into flow cytometry tubes. Wash once with pre-chilled PBS, centrifuge at 300 × g for 5 minutes, and discard the supernatant.
 - 4 Resuspend cells in 100 μL flow cytometry staining buffer for cell counting.
- (2) Surface Staining
 - ① Dosage example for a 100 μL reaction system: CD3-PerCPCy5.5 (1:100), CD56-PE (1:100).
 - ② Add antibodies to the cell tube, gently mix by pipetting, and incubate at 4 °C in the dark for 30 minutes.
 - 3 Wash with 2 mL staining buffer, centrifuge at 300 × g for 5 minutes, and discard supernatant. Repeat washing once.
- (3) Fixation and Membrane Permeabilization
 - ① Resuspend cells in 100 μL fixative, incubate at room temperature in the dark for 15-20 minutes.
 - 2 Add 2 mL staining buffer, centrifuge, discard supernatant, and thoroughly remove fixative.
 - 3 Resuspend cells in 100 μL lysis buffer, incubate at room temperature in the dark for 15-20 minutes.
- (4) Intracellular Staining
 - ① Dilute intracellular antibodies in lysis buffer: Perforin-FITC (1:20), Granzyme-BV421 (1:20). Centrifuge and discard lysis buffer.
 - 2 Add the intracellular antibody cocktail to the cell tube, gently pipette to mix, and incubate at room temperature in the dark for 30-45 minutes.
 - 3 Wash with 2 mL lysis buffer, centrifuge at 300 \times g for 5 minutes, and discard the supernatant. Repeat the wash once.
 - 4 Resuspend cells in 300-500 μL flow cytometry staining buffer and analyze immediately.
- (5) Data Analysis:
 - 1) The suspension was transferred to flow cytometry sample tubes.
 - 2 Analysis strategy:

Gate to exclude doublets (using FSC-A vs FSC-H).

Gate to exclude dead cells (using a viability dye).

Identify the target cell population (NK cells: CD3⁻ CD56⁺).

Analyze PFN and GrB expression: percentage of positive cells (%) and MFI.

Time: 2025.09.10-2025.09.15

6. Result:




Figure.1 PFN, GrB detection of NK cells pre-activated by IL-2, B51G35R, B51G35R-G, B34G35R, B34G35R-G, B51G9-G separately at both pH=6.4 and pH=7.4. Data are representative of at least three independent experiments (**P < 0.01, *P < 0.05).