

Time: 2024.05.10-2024.05.11

1. **Experiment:** Infection of NK92 cell line
2. **Time:** 2024.05.10-2024.05.11
3. **Member:** Hanyue Liu, Xudong Tang, Xiaoyuan Chen
4. **Material:** 5% FBS, complete medium (RPMI-1640 with 10% FBS), polybrene.
5. **Method:**
 - (1) Cell Inoculation:
 - ① Collect NK92 cells from a culture flask.
 - ② Centrifuge at $300\times g$ for 5 minutes at room temperature.
 - ③ Carefully remove the supernatant and resuspend the cell pellet in an appropriate volume of complete medium.
 - ④ Calculate the required volume of cell suspension to achieve a final concentration of 1×10^7 cells/mL.
 - ⑤ Add the calculated volume of cells to an appropriate volume of complete medium to achieve a final concentration of 1×10^7 cells/mL.
 - ⑥ Dispense 100 μL of the 1×10^7 cells/mL suspension into each well of a 96-well plate. Incubate the plates at 37°C in a 5% CO_2 incubator for 24 hrs to achieve approximately 50% confluence using a light microscope to observe the cells. At 50% confluence, the cells should cover approximately half of the well surface. You'll see some areas with dense cell layers and others with less or no cells.
 - (2) Infection:

To each well, 4 μL of a $25\times$ polybrene infection-enhancing solution was added. Considering the multiplicity of infection (MOI) and the virus titer, the specific volume of the virus, calculated as $(\text{MOI} \times \text{cell number}) / \text{virus titer}$, was added to the wells. The plates were then incubated at 37°C for 12-16 hrs.