SynthImmunol_NMU

Time: 2024.06.21-2024.06.22

- 1. Experiment: Infection of CAR on NK92 cell line
- **2.** Time: 2024.06.21-2024.06.22
- 3. Member: Hanyue Liu, Yaqi Gao
- 4. Material: Centrifuge 5810 R (Eppendorf, 5810R), Eclipse TS100 (Nikon, TS100), 5% FBS, complete medium (RPMI-1640 with 10% FBS), polybrene

5. Method:

- (1) Cell Inoculation:
 - ① Collect NK92 cells from a culture flask.
 - (2) Centrifuge at $300 \times g$ for 5 mins at room temperature.
 - ③ Carefully remove the supernatant and resuspend the cell pellet in an appropriate volume of complete medium (e.g., RPMI-1640 with 10% FBS).
 - (4) Calculate the required volume of cell suspension to achieve a final concentration of 1×10^7 cells/mL.
 - (5) Add the calculated volume of cells to an appropriate volume of complete medium to achieve a final concentration of 1×10^7 cells/mL.
 - (6) Dispense 100 μ L of the 1×10⁷ cells/mL suspension into each well of a 96-well plate. Incubate the plates at 37°C in a 5% CO₂ incubator for 24 hrs to achieve approximately 50% confluence using a light microscope to observe the cells. At 50% confluence, the cells should cover approximately half of the well surface. You'll see some areas with dense cell layers and others with less or no cells.
- (2) Infection: To each well, 4 μL of a 25× polybrene infection-enhancing solution was added. Considering the multiplicity of infection (MOI) and the virus titer, the specific volume of the virus, calculated as (MOI × cell number) / virus titer, was added to the wells. The plates were then incubated at 37°C for 12-16 hrs.