Time: 2024.10.26-2024.12.01

- 1. Experiment: Lentivirus Packaging and Infection of CAR on NK-92 cell line
- **2. Time:** 2024.10.26-2024.12.01
- 3. Member: Fan Yang, Qiwen Jiang, Xinxin Zhang, Kaiqing Zhang, Meng Sun

4. Principle:

Lentiviral packaging is a technique to generate replication-incompetent recombinant lentiviral particles for efficient delivery of exogenous genes into target cells. Its core principle relies on the split of viral functional components into separate plasmids, which are co-expressed in packaging cells to assemble infectious but non-replicating viral particles.

5. Materials: Lenti-293 T cells, 5% CO₂ incubator, 5% FBS, CPT Transfection Kit (Viraltherapy), Centrifuge 5810 R (Eppendorf, 5810R), Eclipse TS100 (Nikon, TS100), 96-well plates, 5% CO₂ incubator, complete medium (RPMI-1640 with 10% FBS), polybrene.

6. Method:

- (1) One day before transfection, Lenti-293 T cells was passaged and evenly plated on 100 mm culture dishes (at 8-10×10⁶ cells per dish). Cells were cultured at 37 °C in a 5% CO₂ incubator. The medium was replaced with 10.5 mL of pre-warmed, resistance-free medium (5% FBS) 2-4 hours before transfection.
- (2) A mixture of twenty-three plasmids in equimolar amounts was prepared, and transfection was performed using the calcium phosphate transfection method with the CPT Transfection Kit (Viraltherapy). The transfection mixture composition was as follows:

	Table.1 Transfection mixture composition		
Transfection tube	Reagent	For 1×15cm Dish	For 10×15cm Dish
A	Buffer A	500 μL	45 mL
	lentiviral plasmids	21 μg	210 μg
	pMD2G	7 μg	70 μg
В	pCMVR8.74	14 μL	140 μL
	sterile water	Up to $450~\mu L$	Up to 4.5 mL
	Buffer B	50 μL	500μL

Table.1 Transfection mixture composition

- (3) The homogenous mixtures in tube B were added dropwise into A using a pipette gun and placed for 30 mins at room temperature.
- (4) Cell Inoculation:
 - 1) Q Collect NK-92 cells from a culture flask.
 - 2 Centrifuge at $300 \times g$ for 5 mins at room temperature.
 - ③ Carefully remove the supernatant and resuspend the cell pellet in an appropriate volume of complete medium (e.g., RPMI-1640 with 10% FBS).
 - 4 Calculate the required volume of cell suspension to achieve a final concentration of 1×107 cells/mL.
 - \bigcirc Add the calculated volume of cells to an appropriate volume of complete medium to achieve a final concentration of 1×107 cells/mL.
 - $^{\circ}$ Dispense 100 μL of the 1×107 cells/mL suspension into each well of a 96-well plate. Incubate the plates at 37 $^{\circ}$ C in a 5% CO2 incubator for 24 hours to achieve approximately 50% confluence using a light microscope to observe the cells. At 50% confluence, the cells should cover approximately half of the well surface. You'll see some areas with dense cell layers and others with less or no cells.
- (5) Infection: To each well, 4 μL of a 25 × polybrene infection-enhancing solution was added. Considering

Time: 2024.10.26-2024.12.01

the multiplicity of infection (MOI) and the virus titer, the specific volume of the virus, calculated as $(MOI \times cell \ number) / virus \ titer$, was added to the wells. The plates were then incubated at 37 °C for 12-16 hours.